# SUTURE

# FUTURE

# PLUS

VICRYL PLUS  
MONOCRYL PLUS  
PDS PLUS  
STRATAFIX SPIRAL MONOCRYL PLUS  
STRATAFIX SPIRAL PDS PLUS  
STRATAFIX SYMMETRIC PDS PLUS

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### WHO

The World Health Organization, 2016:

»The panel suggests the use of triclosan-coated sutures for the purpose of reducing the risk of SSI, independent of the type of surgery.« (Conditional recommendation, moderate quality of evidence).<sup>18</sup>

### CDC

Centers for Disease Control and Prevention, 2017:

Guidelines for the Prevention of SSI 2017: »Consider the use of triclosan-coated sutures for the prevention of SSI.«<sup>19</sup>

### EUNETHTA

European Network for Health Technology Assessment, 2017:

»A statistically significant benefit of triclosan-coated sutures in reducing the risk of total incisional SSIs was demonstrated in our SR/MA, based on moderate quality RCTs data«<sup>20</sup>

### ACS & SIS

American College of Surgeons & Surgical Infection Society, 2016:

SSI-Guidelines: »Triclosan antibacterial suture use is recommended for wound closure in clean and clean-contaminated abdominal cases ...«<sup>21</sup>

### NICE

National Institute for Health and Care Excellence, 2013:

»Antimicrobial-coated sutures may reduce surgical site infection risk versus uncoated sutures, although this effect may be specific to particular types of surgery (such as abdominal procedures).«<sup>22</sup>

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**16** De Jonge kommt in seiner Metaanalyse von 2017 (De Jonge, S. W., Atema, J.J., Solomkin, J. S., Boormeester, M. A.. (2017). Meta-analysis and trial sequential analysis of triclosan-coated sutures for the prevention of surgical-site infection. *British Journal of Surgery*. *BJS* 2017; 104: e118–e133. Published online 23 December 2016 in Wiley Online Library ([www.bjs.co.uk](http://www.bjs.co.uk)). DOI: 10.1002/bjs.10445) zu folgenden Ergebnis: »Conclusion: GRADE assessment shows moderate-quality evidence that TCS are effective in reducing SSI. Trial sequential analysis indicates that the effect was robust, and additional data are unlikely to alter the summary effect.« »With the presented moderate quality of evidence based on a large number of trials, it is not thought justified to expose thousands of patients to randomization of an intervention that has been tested so rigorously«. **17** Kommission für Krankenhaushygiene und Infektionsprävention (KRINKO), 2018: Prävention postoperativer Wundinfektionen; S. 461ff, [www.rki.de/DE/Content/Infekt/Krankenhaushygiene/Kommission/Downloads/Empf\\_postopWI.pdf?\\_\\_blob=publicationFile](http://www.rki.de/DE/Content/Infekt/Krankenhaushygiene/Kommission/Downloads/Empf_postopWI.pdf?__blob=publicationFile); Zugriff am 9. 4. 2018. **18** World Health Organization. Guideliness for the prevention of Surgical Site Infection. Geneva, Switzerland, [www.who.int/infection-prevention/publications/ssi-prevention-Guidelines/en](http://www.who.int/infection-prevention/publications/ssi-prevention-Guidelines/en). **19** Berrios-Torres S., Umscheid C., Bratzler D., et al. Centers for Disease Control and Prevention Guidelines for the Prevention of Surgical Site Infection, 2017. *JAMA Surg*. doi:10.1001/jamasurg.2017.0904. **20** [www.eunetha.eu/wp-content/uploads/2018/01/OTCA02\\_Antibacterial-coated-sutures-for-the-prevention-of-abdominal-SSI\\_0.pdf](http://www.eunetha.eu/wp-content/uploads/2018/01/OTCA02_Antibacterial-coated-sutures-for-the-prevention-of-abdominal-SSI_0.pdf) **21** Ban K., Minei J., Laronga C., et al. American College of Surgeons and Surgical Infection Society: Surgical Site Infection Guideliness, 2016 Update. *J Am Coll Surg*. 2016;224:59–74. **22** [www.nice.org.uk/guidance/qs49](http://www.nice.org.uk/guidance/qs49)

## Studien

**1** Ford HR et al. Intraoperative handling and wound healing: controlled clinical trial comparing coated VICRYL Plus antibacterial suture (coated polyglactin 910 suture with triclosan) with coated VICRYL suture (coated polyglactin 910 suture). *Surg Infect (Larchmt)* 2005; 6: 313-21 **2** DeFazio A et al. Does the Use of Vicryl Plus Antibacterial Suture Decrease the Incidence of Umbilical Infection when Compared to Vicryl Suture? *J Minim Invasive Gynecol* 2005; 12 (Suppl.): 37-8 **3** Fleck T et al. Triclosan-coated sutures for the reduction of sternal wound infections: economic considerations. *Ann Thorac Surg* 2007; 84: 232-6 **4** Rozzelle CJ et al. Antimicrobial suture wound closure for cerebrospinal fluid shunt surgery: a prospective, double-blinded, randomized controlled trial. *J Neurosurg Pediatr* 2008; 2: 111-7 **5** Picó RB et al. Prospective study comparing the incidence of wound infection following appendectomy for acute appendicitis in children: conventional treatment versus using reabsorbable antibacterial suture or gentamicin-impregnated collagen fleeces. *Cir Pediatr* 2008; 21: 199-202 **6** Justinger C et al. Antibacterial coating of abdominal closure sutures and wound infection. *Surgery* 2009; 145: 330-4 **7** Zhuang CP et al. Comparison of two absorbable sutures in abdominal wall incision. *Journal of Clinical Rehabilitative Tissue Engineering Research* 2009; 13: 4045-48 **8** Mingmalairak C et al. Efficacy of antimicrobial coating suture coated polyglactin 910 with triclosan (Vicryl Plus) compared with polyglactin 910 (Vicryl) in reduced surgical site infection of appendicitis, double blind randomized control trial, preliminary safety report. *J Med Assoc Thai* 2009; 92: 770-5 **9** Singh H et al. Antibacterial suture reduces surgical site infections in coronary artery bypass grafting. *Heart Surgery Forum* 2010; 13: S85 **10** Zhang ZT et al. Cosmetic outcome and surgical site infection rates of antibacterial absorbable (Polyglactin 910) suture compared to Chinese silk suture in breast cancer surgery: a randomized pilot research. *Chin Med J (Engl)* 2011; 124: 719-24 **11** Chen SY et al. Do antibacterial-coated sutures reduce wound infection in head and neck cancer reconstruction? *Eur J Surg Oncol* 2011; 37: 300-4 **12** Rasić Z et al. Efficacy of antimicrobial triclosan-coated polyglactin 910 (Vicryl\* Plus) suture for closure of the abdominal wall after colorectal surgery. *Coll Antropol* 2011; 35: 439-43 **13** Galal I & El-Hindawy K. Impact of using triclosan-antibacterial sutures on incidence of surgical site infection. *Am J Surg* 2011; 202: 133-8 **14** Justinger C et al. Triclosan-coated sutures

reduce wound infections after hepatobiliary surgery – a prospective nonrandomized clinical pathway driven study. *Langenbecks Arch Surg* 2011; 396: 845-50

**15** Stadler S & Fleck T. Triclosan-coated sutures for the reduction of sternal wound infections? A retrospective observational analysis. *Interact Cardiovasc Thorac Surg* 2011; 13: 296-9

**16** Baracs J et al. Surgical site infections after abdominal closure in colorectal surgery using triclosan-coated absorbable suture (PDS Plus) vs. uncoated sutures (PDS II): a randomized multicenter study. *Surg Infect (Larchmt)* 2011; 12: 483-9

**17** Khachatryan N et al. Prevention of postoperative infections in abdominal surgery using reabsorbable suture with antibacterial activity (Vicryl Plus) versus reabsorbable standard sutures. *Surg Infect (Larchmt)* 2011; 12: A13-A14

**18** Mattavelli I et al. Triclosan-coated suture to reduce surgical site infection after colorectal surgery. *Surg Infect (Larchmt)* 2011; 12: A14-A15

**19** Williams N et al. Randomized trial of antimicrobial-coated sutures to prevent surgical site infection after breast cancer surgery. *Surg Infect (Larchmt)* 2011; 12: 469-74

**20** Isik I et al. Efficiency of antibacterial suture material in cardiac surgery: a double-blind randomized prospective study. *Heart Surg Forum* 2012; 15: E40-5

**21** Seim BE et al. Triclosan-coated sutures do not reduce leg wound infections after coronary artery bypass grafting. *Interact Cardiovasc Thorac Surg* 2012; 15: 411-5

**22** Turtiainen J et al. Effect of triclosan-coated sutures on the incidence of surgical wound infection after lower limb revascularization surgery: a randomized controlled trial. *World J Surg* 2012; 36: 2528-34

**23** Laas E et al. Antibacterial-coated suture in reducing surgical site infection in breast surgery: a prospective study. *Int J Breast Cancer* 2012; 2012: 819578. doi: 10.1155/2012/819578

**24** Hoshino S et al. A study of the efficacy of antibacterial sutures for surgical site infection: a retrospective controlled trial. *Int Surg* 2013; 98: 129-32

**25** Nakamura T et al. Triclosan-coated sutures reduce the incidence of wound infections and the costs after colorectal surgery: a randomized controlled trial. *Surgery* 2013; 153: 576-83

**26** Ueno M et al. Triclosan-coated sutures reduce wound infections after spinal surgery: a retrospective, nonrandomized, clinical study. *Spine J* 2015; 15: 933-8

**27** Justinger C et al. Surgical-site infection after abdominal wall closure with triclosan-impregnated polydioxanone sutures: Results of a randomized clinical pathway facilitated trial

(NCT00998907). *Surgery* 2013; 154: 589-95

**28** Thimour-Bergström L et al. Triclosan-coated sutures reduce surgical site infection after open vein harvesting in coronary artery bypass grafting patients: a randomized controlled trial. *Eur J Cardiothorac Surg* 2013; 44: 931-8

**29** Room H et al. Antibacterial coated sutures reduce laparoscopic postoperative surgical site infections. *Int J Surg* 2013; 11: 722 (Abstract 1330)

**30** Fracalvieri et al. Predictors of wound infection in elective colorectal surgery. Multicenter observational case-control study. *Cir Esp* 2014; 92: 478-84

**31** Jensen CD et al. A Randomized Controlled Trial of Triclosan-Coated Sutures in 2,547 Lower Limb Arthroplasty Operations. *AAOS Annual Meeting* 2014, New Orleans, Louisiana

**32** Diener M et al. Effectiveness of triclosan-coated PDS Plus versus uncoated PDS II sutures for prevention of surgical site infection after abdominal wall closure: the randomised controlled PROUD trial. *Lancet*. 2014; 384: 142-52.

**33** Mattavelli I. et al., Multi-center randomized controlled trial on the effect of triclosan-coated sutures on surgical site infection after colorectal surgery. *Surg Infect (Larchmt)* 2015; 16: 226–235.

**34** Ruiz-Tovar J. et al., Association between triclosan-coated sutures for abdominal wall closure and incisional surgical site infection after open surgery in patients presenting with fecal peritonitis: a randomized clinical trial, *Surgical Infections*. 2016;16(5):588-594.

**35** Steingrimsdóttir S., et al., Triclosan-coated sutures and sternal wound infections: a prospective randomized clinical trial, *European Journal of Clinical Microbiology and Infectious Diseases*. 2015;34(12):2331-2338.

**36** Renko M. et al., Triclosan-containing sutures versus ordinary sutures for reducing surgical site infections in children: a double-blind, randomised controlled trial, *Lancet Infect Dis* 2017; 17: 50–57.

**37a, b** Ruiz-Tovar J. et al., Effect of the Abdominal Fascial Closure with Triclosan-Coated Sutures in Fecal Peritonitis, on Surgical Site Infection, and Evisceration: A retrospective Multi-Center Study, *Surg Infect* 2018; 19: 61-64

## Metaanalysen

**Chang WK**, Srinivasa S, Morton R, et al. Triclosan-Impregnated Sutures to Decrease Surgical Site Infections: Systematic Review and Meta-Analysis of Randomized Trials. *Ann Surg* 2012; 255:854859.

**Sajid MS**, Craciunas L, Sains P, et al. Use of antibacterial sutures for skin closure in controlling surgical site infections: a systematic review of published randomized, controlled trials. *Gastroenterol Rep* 2013;1:4250.

**Edmiston CE**, Daoud FC, Leaper D. Is there an evidence-based argument for embracing an antimicrobial (triclosan)-coated suture technology to reduce the risk for surgical-site infections?: A meta-analysis. *Surgery* 2013;154:89100.

**Wang ZX**, Jiang CP, Cao Y, et al. Systematic review and meta-analysis of triclosan-coated sutures for the prevention of surgical-site infection. *Brit J Surg* 2013;100:465473.

**Daoud FC**. Systematic Literature Review Update of the PROUD Trial: Potential Usefulness of a Collaborative Database. Letter to *Surg Infect (Larchmt)* 2014;15:857858.

**Apisarnthanarak A**, Singh N, Bandong AN, et al. Triclosan-Coated Sutures Reduce the Risk of Surgical Site Infections: A Systematic Review and Meta-Analysis. *Infection Control & Hospital Epidemiology* 2015;36:169179.

**Guo J**, Pan LH, Li YX, et al. Efficacy of triclosan-coated sutures for reducing risk of surgical site infection in adults: a meta-analysis of randomized clinical trials. *J Surg Res* 2016; 201:105117.

**Wu X**, Kubilay NZ, Ren J, et al. Antimicrobial-coated sutures to decrease surgical site infections: a systematic review and meta-analysis. *Eur J Clin Microbiol Infect Dis*. 2016. DOI: 10.1007/s10096-016-2765-y.

**Sandini M**, Mattavelli I, Nespoli L, Uggeri F, Gianotti L. Systematic review and meta-analysis of sutures coated with triclosan for the prevention of surgical site infection after elective colorectal surgery according to the PRISMA statement. *Medicine*. 2016;95:35(e4057).

**De Jonge SW**, Atema JJ, Solomkin JS, Boermeester MA. Meta-analysis and trial sequential analysis of triclosan-coated sutures for the prevention of surgical site infection. *Brit J Surg*. 2017;ePub-DOI: 10.1002/bjs.10445.

**Leaper D.**, Edmiston Jr., C., Holy, C., Meta-analysis of the potential economic impact following introduction of absorbable antimicrobial sutures, *BJS* 2017; 104: e134–e144

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Johnson & Johnson Medical GmbH
Johnson & Johnson Medical Products GmbH
Johnson & Johnson AG Ethicon